

EXHIBIT A



LAL Division PRODUCT SPECIFICATIONS



For PYROSTAR™ reagents, a 0.015 sensitivity is equivalent to a KTA quantitative range of 0.001 - 10. A sensitivity of 0.03 - 0.25 is equivalent to a KTA quantitative range of 0.01 - 10.

1 **PYROSTAR™ ES-F MULTI KIT (2ML)**

4 X 2ML PYROSTAR™ ES-F MULTI plus
1 X CONTROL STANDARD ENDOTOXIN (500ng/vial)

Product will be provided in the format noted above and the following sensitivities: 0.015EU/ml, 0.03EU/ml, 0.06EU/ml, 0.125EU/ml and 0.25EU/ml

2 **PYROSTAR™ ES-F MULTI KIT (5.2ML)**

4 X 5.2ML PYROSTAR™ ES-F MULTI plus
1 X CONTROL STANDARD ENDOTOXIN (500ng/vial)

Product will be provided in the format noted above and the following sensitivities: 0.015EU/ml, 0.03EU/ml, 0.06EU/ml, 0.125EU/ml and 0.25EU/ml

3 **PYROSTAR™ ES-F MULTI (2ML)**

100 X 2ML PYROSTAR™ ES-F MULTI

Product will be provided in the format noted above and the following sensitivities: 0.015EU/ml, 0.03EU/ml, 0.06EU/ml, 0.125EU/ml and 0.25EU/ml

4 **PYROSTAR™ ES-F MULTI (5.2ML)**

100 X 5.2ML PYROSTAR™ ES-F MULTI

Product will be provided in the format noted above and the following sensitivities: 0.015EU/ml, 0.03EU/ml, 0.06EU/ml, 0.125EU/ml and 0.25EU/ml

5 **CONTROL STANDARD ENDOTOXIN**

CONTROL STANDARD ENDOTOXIN 500NG/VIAL

Product will ONLY be provided as sets of 6 x 500ng/vial.

6 **LYSATE REAGENT WATER**

ENDOTOXIN-FREE WATER

Product will ONLY be sold in cases of 12 x 100ml vials or 20 x 30ml vials

7 **GEL CLOT REACTION TUBES**

ENDOTOXIN-FREE GEL CLOT REACTION TUBES

Product will ONLY be provided in cases of 200 tubes; each tube is 10 x 75ml or with aluminum caps in boxes of 50 tubes

